

## Product datasheet for **AP31795HR-N**

### Chicken IgG / Chicken IgY (H+L chain) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA, Ie, WB
Recommended Dilution:	<b>Sandwich ELISA:</b> 1/120,000 dilution of antibody was found to generate an O.D of 1.0 in a 15 minute reaction with Tetramethyl Benzidine as the substrate. <b>Western blot</b> (1/5000). <b>Immunoelectrophoresis.</b>
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified Chicken IgY (whole protein) emulsified in Freund's adjuvant. After multiple injections, goats were bled and serum collected.
Isotype:	Ig
Formulation:	50mM HEPES, pH 7.0, 0.1M NaCl containing 0.2% BSA as stabilizer and 0.01% Thimerosal as preservative Label: HRP State: Liquid purified IgY fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C in the dark.



[View online »](#)