

Product datasheet for **AP31799PU-N**

Rabbit IgG (H+L chain) Chicken Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG (H+L chain) Chicken Polyclonal Antibody
Applications:	ELISA, ID, IF, IHC, WB
Recommended Dilution:	ELISA. Western Blot (1/1000). Immunocytochemistry. Immunohistochemistry (1/250). Quality Control: Antibodies were tested for specificity using Double-Immunodiffusion assays in agarose.
Reactivity:	Rabbit
Host:	Chicken
Immunogen:	Purified, non-denatured Rabbit IgG (containing both heavy and light chains). Production. Eggs were collected from hyperimmunized hens and the IgY fraction purified from the yolks. Anti-Rabbit IgG's were affinity-purified over a column conjugated with purified Rabbit IgG (containing both heavy and light chains).
Isotype:	Ig
Formulation:	PBS, pH 7.2 containing 0.02% Sodium Azide as preservative. State: Aff - Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C in the dark.



[View online »](#)