

Product datasheet for **AP31910PU-N**

MC1 Receptor (MC1R) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Reactivity:	Human, Monkey
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Purified
Concentration:	lot specific
Conjugation:	Unconjugated
Gene Name:	melanocortin 1 receptor
Database Link:	Entrez Gene 4157 Human Q01726
Synonyms:	Melanotropin receptor, Melanocortin receptor 1, MC1R, MSHR, MSH-R



[View online »](#)