

## Product datasheet for **AP32052HR-N**

### Human Lambda Light Chain Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA:</b> 1/4000-1/8000. <b>Western Blot.</b> <b>Immunohistochemistry on Frozen Sections.</b>
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Human lambda light chains.
Specificity:	This antibody recognises Human Lambda Light Chains. It has been cross absorbed with pooled human myeloma proteins with kappa light chains.
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 50% Glycerol Preservative: 0.01% Thiomersal Label: Horseradish Peroxidase
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)