

Product datasheet for **AP32071PU-N**


PDLIM2 Goat Polyclonal Antibody

Product data:

| | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, WB |
| Recommended Dilution: | Peptide ELISA: 1/128000 (Detection Limit). Western blot: 0.3-1 µg/ml. Approx 38kDa bands observed in Mouse Lung lysates (calculated MW of 37.7kDa according to Mouse NP_666090.1). |
| Reactivity: | Bovine, Human, Mouse, Porcine, Rat |
| Host: | Goat |
| Clonality: | Polyclonal |
| Immunogen: | Peptide with sequence from the internal region of the protein sequence according to NP_789847.1; NP_067643.3; NP_932159.1. |
| Specificity: | This antibody is expected to recognize reported isoforms 1 and 2 of PDLIM2 (NP_789847.1; NP_067643.3). |
| Formulation: | Tris saline, pH~7.3 State: Aff - Purified State: Liquid purified Ig fraction Stabilizer: 0.5% BSA Preservative: 0.02% Sodium Azide |
| Concentration: | lot specific |
| Purification: | Ammonium Sulphate Precipitation followed by Antigen Affinity Chromatography using the immunizing peptide |
| Conjugation: | Unconjugated |
| Storage: | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| Gene Name: | PDZ and LIM domain 2 |
| Database Link: | Entrez Gene 64236 Human Q96JY6 |
| Synonyms: | PP6345, PDZ-LIM protein mystique |



[View online »](#)

Product images:

Staining of Mouse Lung lysate using PDLIM2 Antibody at 0.3 ug/ml (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence