

## Product datasheet for **AP32074PU-N**

### Mouse IgG3 Goat Polyclonal Antibody

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Secondary Antibodies   |
| Product Name:         | Mouse IgG3 Goat Polyclonal Antibody  |
| Applications:         | ELISA  |
| Recommended Dilution: | <b>ELISA.</b>  |
| Reactivity:           | Mouse  |
| Host:                 | Goat   |
| Immunogen:            | Mouse IgG3 paraproteins.   |
| Isotype:              | IgG  |
| Formulation:          | Borate buffered saline   |
|                       | State: Purified  |
|                       | State: Liquid purified IgG fraction  |
|                       | Preservative: 0.09% Sodium Azide   |
| Concentration:        | lot specific   |
| Purification:         | Affinity Chromatography  |
| Conjugation:          | Unconjugated   |
| Storage:              | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing. |



[View online »](#)