

Product datasheet for AP32277AP-N

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

Applications: ELISA, IHC, WB

Recommended Dilution: ELISA: 1/2,000-1/4,000.

Immunoblotting.

Immunohistochemistry.

Reactivity: Mouse **Host:** Goat

Immunogen: Pooled antisera from goats hyperimmunized with Mouse IgA paraproteins.

Isotype: lgG

Formulation: 50mM Tris/1mM MgCl2/50% Glycerol, pH 8.0

Label: AP

State: Liquid purifed Ig fraction Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase

Conjugation: AP

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Note: Characterization: To insure lot-to-lot consistency, each batch of product is tested by

ELISA, FLISA, and/or flow cytometry for conformance with characteristics of a standard

reference reagent.

