

Product datasheet for **AP32277HR-N**

Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/4,000-1/8,000. Immunoblotting. Immunohistochemistry.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Pooled antisera from goats hyperimmunized with Mouse IgA paraproteins.
Isotype:	IgG
Formulation:	50% Glycerol/50% PBS, pH 7.4 Label: HRP State: Liquid purified Ig fraction Preservative: None Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Characterization: To insure lot-to-lot consistency, each batch of product is tested by ELISA, FLISA, and/or flow cytometry for conformance with characteristics of a standard reference reagent.



[View online »](#)