

## Product datasheet for **AP32277RP-N**

### Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Goat Polyclonal Antibody
Applications:	ELISA, FC, IF, IHC
Recommended Dilution:	<b>FLISA.</b> <b>Flow Cytometry.</b> <b>Immunofluorescence:</b> ≤ 0.1 µg/10 <sup>6</sup> cells. <b>Immunohistochemistry.</b>
Reactivity:	Mouse
Host:	Goat
Immunogen:	Pooled antisera from goats hyperimmunized with Mouse IgA paraproteins.
Isotype:	IgG
Formulation:	PBS containing 0.09% Sodium Azide and a stabilizing agent Label: PE State: Liquid purified Ig fraction Label: R-Phycoerythrin (R-PE)
Concentration:	lot specific
Conjugation:	PE
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Note:	<b>Characterization:</b> To insure lot-to-lot consistency, each batch of product is tested by ELISA, FLISA, and/or flow cytometry for conformance with characteristics of a standard reference reagent.



[View online »](#)