

## Product datasheet for **AP32277TC-N**

### Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Goat Polyclonal Antibody
Applications:	ELISA, FC, IHC
Recommended Dilution:	<b>FLISA:</b> 1/100-1/400. <b>Flow Cytometry.</b> <b>Immunohistochemistry.</b>
Reactivity:	Mouse
Host:	Goat
Immunogen:	Pooled antisera from goats hyperimmunized with Mouse IgA paraproteins.
Isotype:	IgG
Formulation:	PBS Label: TRITC State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide Label: Rhodamine
Concentration:	lot specific
Conjugation:	TRITC
Storage:	Store undiluted at 2-8°C. This products is photosensitive and should be protected from light.
Note:	<b>Characterization:</b> To insure lot-to-lot consistency, each batch of product is tested by ELISA, FLISA, and/or flow cytometry for conformance with characteristics of a standard reference reagent.



[View online »](#)