

Product datasheet for **AP32277TR-N**

Mouse IgA (alpha chain specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Goat Polyclonal Antibody
Applications:	ELISA, FC, IHC
Recommended Dilution:	FLISA: 1/100-1/400. Flow Cytometry. Immunohistochemistry.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Pooled antisera from goats hyperimmunized with Mouse IgA paraproteins.
Isotype:	IgG
Formulation:	PBS Label: Texas Red State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Texas Red
Storage:	Store undiluted at 2-8°C. This products is photosensitive and should be protected from light.
Note:	Characterization: To insure lot-to-lot consistency, each batch of product is tested by ELISA, FLISA, and/or flow cytometry for conformance with characteristics of a standard reference reagent.



[View online »](#)