

Product datasheet for **AP32845PU-S**

Actin Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Quality Control: This anti-Actin preparation was validated by Western blotting on nitrocellulose membranes (at a concentration of 1 µg/ml) using HRP-Conjugated Goat anti-Chicken IgY (1/5000 dilution, <i>Cat.-No</i> AP31795HR-N) as the secondary reagent. <i>Recommended Dilutions:</i> 1/1000-1/2000 for Western blot analysis of Nitrocellulose membranes. Please note that these dilutions are meant to serve as starting points, and that optimal dilutions may vary. Do Not use this product with PVDF, or other nylon-based membranes, because this often produces unacceptably high backgrounds in chemiluminescence studies.
Reactivity:	Rabbit
Host:	Chicken
Isotype:	IgY
Clonality:	Polyclonal
Immunogen:	Chickens were immunized with actin purified from Rabbit muscle. After repeated injections, immune eggs were collected, and the IgY fractions were purified from the yolks. These IgY fractions were then affinity-purified using a peptide column. The final product contains 50% glycerol (v/v), purified anti-Actin IgY fraction (10 mg/ml) spiked with 25 µg/ml of affinity purified anti-Actin IgY, and then filter-sterilized.
Formulation:	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline (PBS, pH 7.2) State: Aff - Purified State: Liquid purified IgY fraction Preservative: 0.02% (w/v) Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at -20°C This product is photosensitive and should be protected from light.



[View online »](#)

Stability: Shelf life: one year from despatch.