

Product datasheet for **AP32848PU-N**

Zopiclone Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 2.5 µg/ml (EIA)
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Zopiclone-BSA/BTG
Specificity:	Recognizes Zopiclone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate Precipitation Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Background:	Zopiclone is a hypnotic agent used in the treatment of insomnia. Known as a one of the "Z-drugs", they were initially thought to be less addictive and/or habit-forming than benzodiazepines. However, this appraisal has shifted somewhat in the last few years as cases of addiction and habituation have been presented.
Synonyms:	Imovane, Zimovane



[View online »](#)