

Product datasheet for **AP32855PU-N**

PIM2 (23-36) Goat Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Peptide ELISA: 1/64000 (Detection Limit). Western blot: 0.5-1 µg/ml. In transfected HEK293 transiently expressing Human PIM2 (myc and DYKDDDDK tagged), a band of approx. 38kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Immunohistochemistry on Paraffin Sections: 5µg/ml. Paraffin embedded Human Heart.
Reactivity:	Bovine, Human, Mouse, Rat
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence from the internal region (near N Terminus) of the protein sequence according to NP_006866.2.
Specificity:	This antibody recognizes PIM2 (23-36).
Formulation:	Tris saline, pH~7.3 State: Aff - Purified State: Liquid purified IgG fraction Stabilizer: 0.5% BSA Preservative: 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Ammonium Sulphate precipitation followed by Antigen Affinity Chromatography using the immunizing peptide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	Pim-2 proto-oncogene, serine/threonine kinase
Database Link:	Entrez Gene 11040 Human Q9P1W9



[View online »](#)

Synonyms: Pim-2, Pim-2h

Product images:



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human PIM2 with C-terminal MYC tag probed with (0.1ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with (0.1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.