

## Product datasheet for AP33082BT-N

## Sheep IgG (H+L chain) Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Secondary Antibodies

**Product Name:** Sheep IgG (H+L chain) Rabbit Polyclonal Antibody

**Applications:** ELISA, IF, IHC, WB

Recommended Dilution: ELISA.

Dot blot.

Immunoblotting.

Immunocytochemistry.

Immunohistochemistry Paraffin Sections.

Can be used in immunocytochemical and immunohistochemical staining to identify and measure IgG, antigen or antibody, at the cellular and subcellular level by staining of appropriately treated cell and tissue substrates, and to demonstrate circulating antibodies in serodiagnostic microbiology and autoimmune diseases; to identify a specific antigen or immune complex using a reference antibody of sheep origin in the middle layer of the indirect test procedure. As a second step an avidin or streptavidin conjugate of the user's choice has to be used. This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining

and interfere with the specific signal. *Recommended Working Dilutions:* 

Histochemistry and Cytochemical Use: 1/100-1/500.

ELISA and comparable non-precipitating antibody-binding assays: 1/1,000-1/10.000.

Reactivity: Sheep
Host: Rabbit

**Immunogen:** Purified normal IgG isolated from pooled sheep serum.

Freund's complete adjuvant is used in the first step of the immunization procedure.

**Isotype:** lgG



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## Sheep IgG (H+L chain) Rabbit Polyclonal Antibody - AP33082BT-N

Formulation: PBS, pH 7.2

Label: Biotin

State: Lyophilized purified hyperimmune IgG fraction

Stabilizer: None

Preservative: None, as it may interfere with the antibody activity.

Label: N-Hydroxysuccinimidobiotin

Molar radio: Biotin/IgG ~4.2

**Reconstitution Method:** Restore with 1 ml sterile distilled water.

Concentration: 10.0 mg/ml

**Purification:** Hyperimmune antisera with strong activity are selected for fractionation by saltprecipitation

and purification of the IgG fraction by DEAE-chromatography.

Conjugation: Biotin

Storage: Store lyophilized at 2-8°C for 6 months or at -20°C long term.

After reconstitution store the antibody undiluted at 2-8°C for one month

or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.

Note: <u>Adsorption:</u> Not required