

Product datasheet for AP33087BT-N

Porcine IgG (H+L chain) Rabbit Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Porcine IgG (H+L chain) Rabbit Polyclonal Antibody

Applications: ELISA, IF, IHC, WB

Recommended Dilution: ELISA.

Dot blot.

Immunoblotting.

Immunocytochemistry.

Immunohistochemistry Paraffin Sections.

Can be used in immunocytochemical and immunohistochemical staining to identify and measure IgG, antigen or antibody, at the cellular and subcellular level by staining of appropriately treated cell and tissue substrates, and to demonstrate circulating antibodies in serodiagnostic microbiology and autoimmune diseases; to identify a specific antigen or immune complex using a reference antibody of swine origin in the middle layer of the indirect test procedure. As a second step an avidin or streptavidin conjugate of the user's choice has to be used. This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining

and interfere with the specific signal. *Recommended Working Dilutions:*

Histochemistry and Cytochemical Use: 1/100-1/500.

ELISA and comparable non-precipitating antibody-binding assays: 1/1,000-1/10,000.

Reactivity: Porcine
Host: Rabbit

Immunogen: Purified normal IgG isolated from pooled Swine serum.

Freund's complete adjuvant is used in the first step of the immunization procedure.

Isotype: lgG



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Porcine IgG (H+L chain) Rabbit Polyclonal Antibody - AP33087BT-N

Formulation: PBS, pH 7.2

Label: Biotin

State: Lyophilized purified hyperimmune IgG fraction

Stabilizer: None

Preservative: None, as it may interfere with the antibody activity.

Label: N-Hydroxysuccinimidobiotin

Molar radio: Biotin/IgG ~4.9

Reconstitution Method: Restore with 1 ml sterile distilled water.

Concentration: 10.0 mg/ml

Purification: Hyperimmune antisera with strong activity are selected for fractionation by saltprecipitation

and purification of the IgG fraction by DEAE-chromatography.

Conjugation: Biotin

Storage: Store lyophilized at 2-8°C for 6 months or at -20°C long term.

After reconstitution store the antibody undiluted at 2-8°C for one month

or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.

Note: <u>Adsorption:</u> Not required.

<u>Conjugation procedure:</u> A proprietary technique for the binding to biotin is used, followed by several purification steps. After each step activity and specificity are tested in a variety of

techniques. The conjugate is lyophilized to assure stability and long shelf life.