

## Product datasheet for **AP33323SU-S**

### Okadaic Acid Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA</b> (1/10,000).
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Okadaic Acid coupled to KLH
Specificity:	Okadaic Acid.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Okadaic acid is a toxin that accumulates in bivalves and causes diarrhetic shellfish poisoning. Okadaic acid was named from the marine sponge <i>Halichondria okadai</i> , from which okadaic acid was isolated for the first time.



[View online »](#)