

Product datasheet for **AP33329SU-N**

Relaxin 2 (RLN2) (1-8) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Relaxin 1B chain (aa 1-9) KLH conjugated
Specificity:	This antibody detects Human Relaxin 1 B chain (aa 1-9). There was no cross reactivity obtained with Human Relaxin 1A chain.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	relaxin 2
Database Link:	Entrez Gene 6019 Human P04090
Synonyms:	Relaxin A chain, Relaxin B chain, RLN1, RLN2, RLN, Prorelaxin H1, Prorelaxin H2



[View online »](#)