

Product datasheet for **AP33332SU-N**

VWA2 (496-732) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry on Paraffin Sections.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant Human AMACO (aa 496-732), C-terminal <i>His-tag</i> .
Specificity:	This antibody detects Recombinant human Amaco (aa 496-732).
Formulation:	State: Serum State: Lyophilised Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	von Willebrand factor A domain containing 2
Database Link:	Entrez Gene 340706 Human Q5GFL6
Synonyms:	CCSP-2



[View online »](#)