

Product datasheet for **AP33333PU-N**

Ne-Carboxymethyllysine (CML) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/200-1/20,000.
Host:	Goat
Clonality:	Polyclonal
Immunogen:	N ^ε -(Carboxymethyl)lysine KLH-conjugated
Specificity:	This antibody specifically binds to CML modified proteins.
Formulation:	75mM PBS, 75mM NaCl, 0.5 mM EDTA, pH 7.2 State: Aff - Purified State: Lyophilized purified Serum Preservative: 0.02% Sodium Azide
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Sepharose Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)