

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for AP33339PU-N

### Amyloid beta A4 protein / APP (556-575) Chicken Polyclonal Antibody

#### **Product data:**


Product Type:	Primary Antibodies
Applications:	IF, IHC
Recommended Dilution:	<ul> <li>Immunocytochemistry.</li> <li>Immunohistochemistry.</li> <li><i>Recommended Dilutions:</i> 1/500-1/1000 for Immunohistochemistry and</li> <li>Immunocytochemistry using 2% paraformaldehyde-fixed tissues or cells.</li> <li><i>Quality Control:</i> This anti-Peptide antibody mixture was analyzed by Immunohistochemistry</li> <li>(1/1000 dilution) using Fluorescein-labeled Goat anti-Chicken IgY (1/500 dilution, CatNo</li> <li>AP31795FC-N) as the secondary reagent.</li> </ul>
Reactivity:	Human, Mouse, Rat
Host:	Chicken
lsotype:	IgY
Clonality:	Polyclonal
Immunogen:	Chickens were immunized with synthetic peptide APP#5 conjugated to keyhole limpet hemocyanin (KLH). The sequence of APP peptide #5 is EEI QDE VDE LLQ KEQ NYS DD, which corresponds to residues #556-575 of the Human gene product (NP_000475.1)]. Identical sequences are also found in the Mouse and Rat homologs of APP.
Specificity:	Recognizes Amyloid Precursor Protein (APP), Peptide #5.
Formulation:	Sodium Phosphate-Buffered (10 mM) isotonic (0.9%, w/v) Saline (PBS, pH 7.2) State: Aff - Purified State: Liquid purified (filter sterilized) IgY fraction Preservative: 0.02% (w/v) Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography using a peptide column. <b>Preparation:</b> After repeated injections, immune eggs were collected, and the IgY fractions were purified from the yolks. These IgY fractions were then affinity-purified using a peptide column, and the concentrations of the eluates adjusted to 100 ug/ml. Finally, equal volumes of these affinity-purified anti-peptide antibodies were mixed, and the preparation was filter- sterilized.
Conjugation:	Unconjugated



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Amyloid beta A4 protein / APP (556-575) Chicken Polyclonal Antibody – AP33339PU-N
Storage:	Store the antibody undiluted in the dark at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	amyloid beta precursor protein
Database Link:	<u>Entrez Gene 351 Human</u> <u>P05067</u>
Synonyms:	Alzheimer disease amyloid protein, Amyloid Precursor Protein, ABPP, APPI, PreA4, Cerebral vascular amyloid peptide, CVAP

## Product images:



Amyloid Precursor Protein (APP), Peptide #5:Mouse cortical neurons in culture. The green staining is APP-immunoreactivity, using Fluorescein-labeled Goat anti-Chicken IgY (Cat.-No [AP31795FC-N]) and Rhodamine-labeled phalloidin as a counterstain. Note the AP

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US