

Product datasheet for **AP33374SU-S**

L-Kynurenic Acid Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/8000.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	L-Kynurenic Acid KLH-conjugated
Specificity:	This antibody detects L-Kynurenic acid. Cross reactivity obtained with 3-Hydroxy-Kynurenine, was less than 0.2%. There was no cross reactivity obtained with Indole Serotonin, or Quinoline.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



Background:	Kynurenic acid (KYNA or KYN) is a product of the normal metabolism of amino acid L-tryptophan. It has been shown that kynurenic acid possesses neuroactive activity. It acts as an antiexcitotoxic and anticonvulsant, most likely through acting as an antagonist at excitatory amino acid receptors. Because of this activity, it may influence important neurophysiological and neuropathological processes. As a result, kynurenic acid has been considered for use in therapy in certain neurobiological disorders. Conversely, increased levels of kynurenic acid have also been linked to certain pathological conditions. Kynurenic acid was discovered in 1853 by the German chemist Justus von Liebig in dog urine, which it was apparently named after. (1) It is formed from L-kynurenine in a reaction catalyzed by the enzyme kynurenine—oxoglutarate transaminase.
Synonyms:	kynuronic acid, quinurenic acid, transtorine, KynA, 4-hydroxyquinoline-2-carboxylic acid, kinurenic acid
Note:	CAS numbe: 492-27-3 Molecular formula: C₁₀H₇NO₃ Molar mass: 189.168 g/mol Melting point: 282.5°C