

Product datasheet for **AP33451FC-N**

Mouse IgG Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG Goat Polyclonal Antibody
Applications:	IF
Recommended Dilution:	Immunofluorescence: 1/500-1/1000. Immunohistochemistry.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Purified Mouse IgG fraction
Isotype:	IgG
Formulation:	PBS, pH 7.4 Label: FITC State: Aff - Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol
Concentration:	lot specific
Purification:	Affinity Chromatorgraphy on Protein A
Conjugation:	FITC
Storage:	Upon receipt, store (in aliquots) at -20°C to -80°C. Avoid repeated freezing and thawing.



[View online »](#)