

Product datasheet for **AP33451HR-N**

Mouse IgG Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/10000. Western Blot: 1/5000-1/10000, Immunohistochemistry.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Purified Mouse IgG fraction
Isotype:	IgG
Formulation:	PBS, pH 7.4 Label: HRP State: Aff - Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	HRP
Storage:	Upon receipt, store (in aliquots) at -20°C to -80°C. Avoid repeated freezing and thawing.



[View online »](#)