

Product datasheet for AP33525PU-M

HAS3 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/2000. Immunohistochemistry on Paraffin Sections: 1/50-1/200.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant full length Human HAS3.
Specificity:	This antibody detects endogenous levels of HAS3 protein.
Formulation:	Phosphate buffered saline (PBS), 50% Glycerol State: Aff - Purified State: Liquid purified lg fraction (> 95% pure by SDS-PAGE). Preservative: 0.09% Sodium Azide
Concentration:	1.0 mg/ml
Purification:	Affinity Chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~63 kDa
Gene Name:	hyaluronan synthase 3
Database Link:	<u>Entrez Gene 3038 Human</u> <u>O00219</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE HAS3 Rabbit Polyclonal Antibody – AP33525PU-M

Background:HAS1, HAS2 and HAS3 are HA (hyaluronan or hyaluronic acid) synthase proteins. The
extracellular matrix in most vertebrates express HA, which is a high molecular weight linear
polysaccharide composed of alternating glucuronic acid and N-acetylglucosamine residues
linked by b-1,3 and b-1,4 glycosidic bonds. The three HAS genes show distinct patterns of
expression during development and their protein products play significantly different roles in
the formation of the HA matrix. Both HAS1 and HAS2 synthesize high molecular-weight HA,
whereas HAS3 produces lower molecular weight HA. The expression of the three HAS
isoforms is more prominent in growing cells than in resting cells and is differentially
regulated by various stimuli, suggesting distinct functional roles of the three proteins. HAS3
produces both secreted and cell-associated forms of hyaluronan and is the most active of the
three isoforms of this enzyme in adults. HAS3 gene expression plays a crucial role in the
regulation of hyaluronan synthesis in the epidermis. Specifically, IFN-g markedly upregulates
HAS3 mRNA, whereas TGFβ downregulates HAS3 transcript levels. The human HAS3 gene
maps to chromosome 16q22.1.

Synonyms: Hyaluronan synthase 3, Hyaluronate synthase 3

Product images:



WesternBlot (WB) analysis of HAS3 polyclonal antibody

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., <u>9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US</u>