

## Product datasheet for AP50072PU-N

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

## **ACTN3 (Center) Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

**Applications:** IHC, WB

Recommended Dilution: ELISA: 1/1000.

Western blotting: 1/100 - 1/500. Immunohistochemistry: 1/10 - 1/50.

Reactivity: Human
Host: Rabbit

**Isotype:** lg

**Clonality:** Polyclonal

Immunogen: KLH conjugated synthetic peptide between 408-437 amino acids from the Central region of

human ACTN3

**Specificity:** This antibody reacts to Alpha-actinin-3.

**Formulation:** PBS containing 0.09% (W/V) sodium azide as preservative

State: Aff - Purified

State: Liquid purified Ig fraction

**Concentration:** lot specific

**Purification:** Affinity chromatography on Protein A

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: actinin alpha 3 (gene/pseudogene)

**Database Link:** Entrez Gene 89 Human

Q08043

Background: Alpha-actinin is an actin-binding protein with multiple roles in different cell types. This gene

expression is limited to skeletal muscle. It is localized to the Z-disc and analogous dense

bodies, where it helps to anchor the myofibrillar actin filaments.

**Synonyms:** Alpha actinin skeletal muscle isoform 3






Note: Molecular Weight: 103241 Da

**Protein Pathways:** Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Focal adhesion,


Leukocyte transendothelial migration, Regulation of actin cytoskeleton, Systemic lupus

erythematosus, Tight junction

## **Product images:**



ACTN3 Antibody (Center) western blot analysis in MDA-MB231 cell line lysates (35ug/lane). This demonstrates the ACTN3 antibody detected the ACTN3 protein (arrow).



ACTN3 Antibody (Center) immunohistochemistry analysis in formalin fixed and paraffin embedded human skeletal muscle followed by peroxidase conjugation of the secondary antibody and DAB staining. This data demonstrates the use of ACTN3 Antibody (Center) for immunohistochemistry. Clinical relevance has not been evaluated.