

Product datasheet for **AP50103PU-N**

AFF4 (N-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western blotting: 1/100 - 1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 1~30 amino acids from the N-terminal region of human AF5q31
Specificity:	This antibody reacts to AFF4.
Formulation:	PBS containing 0.09% (W/V) sodium azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	AF4/FMR2 family member 4
Database Link:	Entrez Gene 27125 Human Q9UHB7
Synonyms:	AF5Q31, HSPC092, AF4/FMR2 family member 4
Note:	Molecular Weight: 127459 Da
Protein Families:	Transcription Factors



[View online »](#)

Product images:



Western blot analysis of AF5q31 Antibody (N-term) in HL-60 cell line lysates (35ug/lane). AF5q31 (arrow) was detected using the purified Pab.