

## Product datasheet for **AP50191PU-N**

### **ANKS4B (Center) Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	<b>ELISA:</b> 1/1000. <b>Western blotting:</b> 1/100 - 1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 183-213 amino acids from the Central region of human ANKS4B
Specificity:	This antibody reacts to ANKS4B.
Formulation:	PBS containing 0.09% (W/V) sodium azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	ankyrin repeat and sterile alpha motif domain containing 4B
Database Link:	<a href="#">Entrez Gene 257629 Human Q8N8V4</a>
Synonyms:	HARP
Note:	<b>Molecular Weight:</b> 51187 Da



[View online »](#)

## Product images:



ANKS4B Antibody (Center) western blot analysis in WiDr cell line lysates (35ug/lane). This demonstrates the ANKS4B antibody detected the ANKS4B protein (arrow).