

Product datasheet for **AP50371PU-N**

Osteocalcin (BGLAP) (C-term)/(C-term C80) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western blotting: 1/1000.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 66-100 amino acids from the C-terminal region of human BGLAP
Specificity:	This antibody recognizes Osteocalcin.
Formulation:	PBS containing 0.09% (W/V) Sodium Azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Saturated Ammonium Sulfate (SAS) precipitation
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	bone gamma-carboxyglutamate protein
Database Link:	Entrez Gene 632 Human P02818
Synonyms:	Bone Gla protein, BGLAP
Note:	Molecular Weight: 10963 Da
Protein Families:	ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)

Product images:



BGLAP Antibody (C80) western blot analysis in HepG2 cell line lysates (35ug/lane). This demonstrates the BGLAP antibody detected the BGLAP protein (arrow).