

Product datasheet for **AP50788PU-N**


CCDC71 (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western blot: 1/100 - 1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 349-379 amino acids from the C-terminal region of human CCDC71
Specificity:	This antibody reacts to CCDC71.
Formulation:	PBS State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.09% (W/V) sodium azide
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	49618 Da
Gene Name:	coiled-coil domain containing 71
Database Link:	Entrez Gene 64925 Human Q8IV32
Synonyms:	Coiled-coil domain-containing protein 71



[View online »](#)

Product images:

CCDC71 Antibody (C-term) western blot analysis in K562 cell line lysates (35ug/lane). This demonstrates the CCDC71 antibody detected the CCDC71 protein (arrow).