

Product datasheet for **AP51279PU-N**


DMRT3 (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western Blot: 1/100-1/500.
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 438-167 amino acids from the C-terminal region of human DMRT3
Specificity:	This antibody recognizes Human and Mouse DMRT3 (C-term).
Formulation:	PBS State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	doublesex and mab-3 related transcription factor 3
Database Link:	Entrez Gene 240590 Mouse Entrez Gene 58524 Human Q9NQL9
Background:	May regulate transcription during sexual development (By similarity)
Synonyms:	Doublesex- and mab-3-related transcription factor 3, DMRTA3
Note:	Molecular Weight: 51199 Da



[View online »](#)

Product images:

Western blot analysis of DMRT3 Antibody (C-term) in mouse kidney tissue lysates (35ug/lane). This demonstrates the DMRT3 antibody detected the DMRT3 protein (arrow).