

Product datasheet for **AP51332PU-N**


DTWD1 (Center) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western Blot: 1/100-1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 151-181 amino acids from the Central region of human DTWD1
Specificity:	This antibody recognizes Human DTWD1 (Center).
Formulation:	PBS containing 0.09% (W/V) Sodium Azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	DTW domain containing 1
Database Link:	Entrez Gene 56986 Human Q8N5C7
Synonyms:	MDS009
Note:	Molecular Weight: 35248 Da



[View online »](#)

Product images:

Western blot analysis in of DTWD1 Antibody (Center) in MDA-MB435 cell line lysates (35ug/lane). This demonstrates the DTWD1 antibody detected the DTWD1 protein (arrow).