

# **Product datasheet for AP51398PU-N**

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

## **EIF3K (C-term) Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

Applications: WB

Recommended Dilution: ELISA: 1/1000.

Western Blot: 1/100 - 1/500.

Reactivity: Human Host: Rabbit

Isotype: lg

Clonality: Polyclonal

Immunogen: KLH conjugated synthetic peptide between 175-204 amino acids from the C-terminal region

of human EIF3K

**Specificity:** This antibody reacts to EIF3K.

Formulation: PBS

State: Aff - Purified

State: Liquid purified Ig fraction

Preservative: 0.09% (W/V) sodium azide

**Concentration:** lot specific

**Purification:** Affinity chromatography on Protein A

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

**Predicted Protein Size:** 25060 Da

**Gene Name:** eukaryotic translation initiation factor 3 subunit K

**Database Link:** Entrez Gene 27335 Human

Q9UBQ5





#### EIF3K (C-term) Rabbit Polyclonal Antibody - AP51398PU-N


Background:

The 700-kD eukaryotic translation initiation factor-3 (eIF3) is the largest eIF and contains at least 12 subunits, including EIF2S12. eIF3 plays an essential role in translation by binding directly to the 40S ribosomal subunit and promoting formation of the 40S preinitiation complex (Mayeur et al., 2003 [PubMed 14519125]).

Synonyms:

Eukaryotic translation initiation factor 3 subunit K, eIF-3 p25, eIF-3 p28, PLAC-24

## **Product images:**



EIF3K Antibody (C-term) western blot analysis in T47D cell line lysates (35ug/lane). This demonstrates the EIF3K antibody detected the EIF3K protein (arrow).