

Product datasheet for **AP51692PU-N**

FNDC7 (C-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	FC, IHC, WB
Recommended Dilution:	ELISA: 1/1000. Western Blot: 1/100-1/500. Flow Cytometry: 1/10-1/50. Immunohistochemistry on Paraffin Sections: 1/10-1/50.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 705-734 amino acids from the C-terminal region of Human FNDC7.
Specificity:	This antibody recognizes Human FNDC7 (C-term).
Formulation:	PBS containing 0.09% (W/V) Sodium Azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A column, followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	fibronectin type III domain containing 7
Database Link:	Entrez Gene 163479 Human Q5VTL7
Synonyms:	FLJ35838; MGC163293; RP11-293A10.2
Note:	Molecular Weight: 78494 Da




[View online »](#)


Product images:



Western blot analysis of FNDC7 Antibody (C-term) in HepG2 cell line lysates (35ug/lane). This demonstrates the FNDC7 antibody detected the FNDC7 protein (arrow).



Immunohistochemistry analysis in formalin fixed and paraffin embedded human rectum tissue reacted with FNDC7 Antibody (C-term) followed by peroxidase conjugation of the secondary antibody and DAB staining.



Flow cytometric analysis of K562 cells using FNDC7 Antibody (C-term) (right histogram) compared to a negative control cell (left histogram). FITC-conjugated donkey-anti-rabbit secondary antibodies were used for the analysis.