

Product datasheet for **AP51744PU-N**


FOXRED1 (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western Blot: 1/100-1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 355~384 amino acids from the C-terminal region of Human FOXRED1
Specificity:	This antibody recognizes Human FOXRED1 (C-term).
Formulation:	PBS containing 0.09% (W/V) Sodium Azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A column, followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	FAD-dependent oxidoreductase domain containing 1
Database Link:	Entrez Gene 55572 Human Q96CU9
Synonyms:	FP634 , H17
Note:	Molecular Weight: 53812 Da



[View online »](#)

Product images:

Western blot analysis of FOXRED1 Antibody (C-term) in NCI-H460 cell line lysates (35ug/lane). FXRD1 (arrow) was detected using the purified Pab.