

Product datasheet for **AP52954PU-N**

NUAK2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western blot: 1/100 - 1/500.
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	Recombinant protein from human NUAK2.
Specificity:	This antibody reacts to NUAK2.
Formulation:	PBS State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.09% (W/V) sodium azide
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	69612 Da
Gene Name:	NUAK family kinase 2
Database Link:	Entrez Gene 81788 Human Q9H093
Synonyms:	SNARK
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)

Product images:



NUAK2 Antibody western blot analysis in mouse NIH-3T3 cell line lysates (35ug/lane). This demonstrates the NUAK2 antibody detected the NUAK2 protein (arrow).



NUAK2 Antibody western blot analysis in NCI-H292 cell line lysates (35ug/lane). This demonstrates the NUAK2 antibody detected the NUAK2 protein (arrow).