

Product datasheet for **AP53914PU-N**


SIAH3 (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1000. Western Blot: 1/100-1/500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 188-217 amino acids from the C-terminal region of Human SIAH3.
Specificity:	This antibody recognizes Human SIAH3 (C-term).
Formulation:	PBS State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	siah E3 ubiquitin protein ligase family member 3
Database Link:	Entrez Gene 283514 Human Q8IW03
Synonyms:	Seven in absentia homolog 3
Note:	Molecular Weight: 30660 Da



[View online »](#)

Product images:

Western blot analysis of SIAH3 (C-term) Antibody in U251 cell line lysates (35ug/lane). This demonstrates the SIAH3 antibody detected the SIAH3 protein (arrow).