

Product datasheet for AP54311PU-N

TNFRSF1A (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: FC, IF, IHC, WB **Recommended Dilution:** ELISA: 1:1:000.

Western blot: 1:100~500.

Immunohistochemistry on paraffin sections: 1:50~100.

Immunofluorescence: 1:10~50. Flow Cytometry: 1:10-1:50.

Reactivity: Human Host: Rabbit

Isotype: lg

Clonality: Polyclonal

Immunogen: KLH conjugated synthetic peptide between 24-52 amino acids from the N-terminal region of

human TNFRSF1A

Specificity: This antibody detects CD120a / TNFR1 (N-term).

Formulation: PBS with 0.09% (W/V) sodium azide

State: Aff - Purified State: Liquid Ig fraction

Concentration: lot specific

Purification: Protein A column followed by peptide affinity purification

Conjugation: Unconjugated

Storage: Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated

freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: tumor necrosis factor receptor superfamily member 1A

Database Link: Entrez Gene 7132 Human

P19438



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:

The protein encoded by this gene is a member of the TNF-receptor superfamily. This protein is one of the major receptors for the tumor necrosis factor-alpha. This receptor can activate NF-kappaB, mediate apoptosis, and function as a regulator of inflammation. Antiapoptotic protein BCL2-associated athanogene 4 (BAG4/SODD) and adaptor proteins TRADD and TRAF2 have been shown to interact with this receptor, and thus play regulatory roles in the signal transduction mediated by the receptor. Germline mutations of the extracellular domains of this receptor were found to be associated with the autosomal dominant periodic fever syndrome. The impaired receptor clearance is thought to be a mechanism of the disease.

Synonyms: Tumor necrosis factor receptor 1, TNF-R1, TNF-R1, TNFR-I, p55, p60, Tnfrsf1a


Note: Molecular Weight: 50495 Da

Protein Families: Druggable Genome, Secreted Protein, Transcription Factors, Transmembrane


Protein Pathways: Adipocytokine signaling pathway, Alzheimer's disease, Amyotrophic lateral sclerosis (ALS),

Apoptosis, Cytokine-cytokine receptor interaction, MAPK signaling pathway

Product images:




TNFRSF1A Antibody (N-term) western blot analysis in U251 cell line lysates (35 ug/lane). This demonstrates the TNFRSF1A antibody detected the TNFRSF1A protein (arrow).



TNFRSF1A Antibody (N-term) immunohistochemistry analysis in formalin fixed and paraffin embedded human rectum tissue followed by peroxidase conjugation of the secondary antibody and DAB staining. This data demonstrates the use of TNFRSF1A Antibody (N-term) for immunohistochemistry. Clinical relevance has not been evaluated.





Confocal immunofluorescent analysis of TNFRSF1A Antibody (N-term) with U-251MG cell followed by Alexa Fluor 488-conjugated goat antirabbit IgG (green). DAPI was used to stain the cell nuclear (blue).