

Product datasheet for **AP54448PU-N**

UBTD1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1;000. Western blot: 1:100~500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	Recombinant protein from human UBTD1
Specificity:	This antibody detects UBTD1.
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	ubiquitin domain containing 1
Database Link:	Entrez Gene 80019 Human Q9HAC8
Synonyms:	Ubiquitin domain-containing protein 1
Note:	Molecular Weight: 25938 Da
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Western blot analysis of UBTD1 (arrow) using rabbit polyclonal UBTD1 Antibody. 293 cell lysates (2 ug/lane) either nontransfected (Lane 1) or transiently transfected (Lane 2) with the UBTD1 gene.