

Product datasheet for **AP54468PU-N**

LOC100128816 (N-term) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1:1,000. Western blot: 1:100~500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 21-50 amino acids from the N-terminal region of human UNQ3104
Specificity:	This antibody detects UNQ3104 (N-term).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff – Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 – 8 °C for up to six months or (in aliquots) at –20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Entrez Gene 100128816 Human Q6UY24
Synonyms:	ACAH3104



Note:

Molecular Weight: 12;518 Da

Product images:



UNQ3104 Antibody (N-term) western blot analysis in MDA-MB231 cell line lysates (35 ug/lane). This demonstrates the UNQ3104 antibody detected the UNQ3104 protein (arrow).