

Product datasheet for **AP54593PU-N**

YDJC (Center) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	ELISA: 1:1,000. Western blot: 1:100~500 Immunohistochemistry on paraffin sections: 1:50~100.
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 176~206 amino acids from the Central region of human YDJC
Specificity:	This antibody detects YDJC (Center).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	YdjC homolog (bacterial)
Database Link:	Entrez Gene 69101 Mouse Entrez Gene 150223 Human A8MPS7




Synonyms: Carbohydrate deacetylase

Note: **Molecular Weight:** 34466 Da

Product images:



YDJC Antibody (Center) western blot analysis in mouse testis tissue lysates (35 ug/lane). This demonstrates the YDJC antibody detected the YDJC protein (arrow).



YDJC Antibody (Center) immunohistochemistry analysis in formalin fixed and paraffin embedded human tonsil tissue followed by peroxidase conjugation of the secondary antibody and DAB staining. This data demonstrates the use of the YDJC Antibody (Center) for immunohistochemistry. Clinical relevance has not been evaluated.