

Product datasheet for **AP54608PU-N**

ZBBX (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	FC, IHC, WB
Recommended Dilution:	ELISA: 1:1,000. Western blot: 1:100~500. Immunohistochemistry: 1:50~100. Flow cytometry: 1:10~50.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 41-70 amino acids from the N-terminal region of human ZBBX
Specificity:	This antibody detects ZBBX (N-term).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	zinc finger B-box domain containing




Database Link: [Entrez Gene 79740 Human A8MT70](#)


Synonyms: FLJ23049

Note: **Molecular Weight:** 91426 Da


Product images:



ZBBX Antibody (N-term) western blot analysis in MDA-MB435 cell line lysates (35 ug/lane). This demonstrates the ZBBX antibody detected the ZBBX protein (arrow).



ZBBX antibody (N-term) immunohistochemistry analysis in formalin fixed and paraffin embedded human colon carcinoma followed by peroxidase conjugation of the secondary antibody and DAB staining. This data demonstrates the use of the ZBBX antibody (N-term) for immunohistochemistry. Clinical relevance has not been evaluated.



ZBBX Antibody (N-term) flow cytometric analysis of MDA-MB435 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis