

Product datasheet for **AP54627PU-N**

RTL3 (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1:1;000. Western blot: 1:100~500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 45-75 amino acids from the N-terminal region of human ZCCHC5
Specificity:	This antibody detects ZCCHC5 (N-term).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	zinc finger CCHC-type containing 5
Database Link:	Entrez Gene 203430 Human Q8N8U3
Synonyms:	FLJ38865; Mar3; Mart3; ZHC5
Note:	Molecular Weight: 52817 Da



[View online »](#)