

Product datasheet for **AP54684PU-N**

ZNF28 (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1:1,000. Western blot: 1:100~500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 87-115 amino acids from the N-terminal region of human ZNF28
Specificity:	This antibody detects ZNF28 (N-term).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff – Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 – 8 °C for up to six months or (in aliquots) at –20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	zinc finger protein 28
Database Link:	Entrez Gene 7576 Human P17035




Background: May be involved in transcriptional regulation.

Synonyms: Zinc finger protein 28, KOX24

Note: **Molecular Weight:** 83658 Da

Product images:



ZNF28 Antibody (N-term) western blot analysis in NCI-H460 cell line lysates (35 ug/lane). This demonstrates the ZNF28 antibody detected the ZNF28 protein (arrow).