

Product datasheet for **AP54729PU-N**


ZNF843 (Center) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1:1;000. Western blot: 1:100~500.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide between 111-141 amino acids from the Central region of human ZNF843
Specificity:	This antibody detects ZNF843 (Center).
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein A column followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to six months or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	zinc finger protein 843
Database Link:	Entrez Gene 283933 Human Q8N446
Synonyms:	Zinc finger protein 843
Note:	Molecular Weight: 37065 Da



[View online »](#)

Product images:


ZNF843 Antibody (Center) western blot analysis in NCI-H292 cell line lysates (35 ug/lane). This demonstrates the ZNF843 antibody detected the ZNF843 protein (arrow).