

Product datasheet for **AP54745SU-N**

14-3-3 sigma (SFN) pSer248 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from C-term domain of the human Stratifin protein.
Specificity:	This antibody reacts with the 27 kDa Stratifin pS248 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	stratifin
Database Link:	Entrez Gene 2810 Human P31947
Synonyms:	Stratifin, HME1



[View online »](#)