

Product datasheet for **AP54756SU-N**

SSH3BP1 (ABI1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from C- term domain of the human Abi-1 protein.
Specificity:	This antibody reacts with 55 kDa Abi-1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	abl interactor 1
Database Link:	Entrez Gene 10006 Human Q8IZP0
Synonyms:	Abl interactor 1, Abelson interactor 1, Abi-1, e3B1, Spectrin binding protein, Nap1-binding protein, Abl-binding protein 4, AblBP4



[View online »](#)