

Product datasheet for **AP54769SU-N**

Aequorin Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/1000-1/10000.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Purified Aequorin protein.
Specificity:	This antibody reacts with 22 kDa Aequorin protein from Jeffyfish.
Formulation:	State: Serum State: Lyophilized Serum Preservative: None
Reconstitution Method:	Restore in 0.2 ml distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P07164



[View online »](#)