

Product datasheet for **AP54780PU-N**

Androgen Receptor (AR) pSer515 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-terminal domain of Human AR (containing phosphorylated Ser515)
Specificity:	This antibody reacts with 99 kDa AR protein with a phosphoserine 515.
Formulation:	0.1M Tris 0.1M Glycine, 2% Sucrose State: Purified State: Lyophilized purified powder Preservative: None
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	androgen receptor
Database Link:	Entrez Gene 367 Human P10275
Synonyms:	Dihydrotestosterone receptor, DHTR, NR3C4



[View online »](#)