

## Product datasheet for **AP54788SU-N**

### Apolipoprotein A V (APOA5) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/100 - 1/500. <b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of the human ApoA5 protein.
Specificity:	This antibody reacts with ~41 kDa Apo-AV.
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	apolipoprotein A5
Database Link:	<a href="#">Entrez Gene 116519 Human Q6Q788</a>
Synonyms:	Apolipoprotein A-V, Apo-AV, ApoA-V, RAP3



[View online »](#)